


# Plan Parameters Report

<b>Hospital Name:</b>	Hospital de Clínicas Caracas
<b>Patient Name:</b>	Abasalo Verdes, Sandra Patricia
<b>Patient ID / ID2 / SSN:</b>	10291238 / 15-396 /
<b>Date of Birth:</b>	-
<b>Course ID / Plan ID:</b>	Curso1 / No_GTV
<b>Plan Created:</b>	Jugleys, martes, 02 de agosto de 2016 11:19 a.m.
<b>Plan Last Modified:</b>	, martes, 02 de agosto de 2016 11:19 a.m.
<b>Plan Approval:</b>	Unapproved, martes, 02 de agosto de 2016 11:19 a.m.
<b>Image ID:</b>	CT_SRS_VMS (image has been processed with Head Frame Detection)
<b>Treatment Unit:</b>	Trilogy
<b>Energy / Dose Rate:</b>	6X-SRS / 1000 MU/min
<b>Fixation:</b>	VMS-Frame
<b>Software Version:</b>	Cone Planning 10.0.42.24719

## Isocenter 1 - Transversal



## Calculation

**Calculation Model:** CDC\_10.0.28.2  
**Calculation Options:** Grid: Sparse  
**Arc Angle Resolution:** 10°  
**Dose Grid Resolution:** 1.0mm  
**Calculation Log:** [SRS]  
 Information: Service: ECDC.SRSDose.10.0.28.9833  
 Information: Servant: p2692@e5620\*  
 Information: Client workstation E5620/3080  
 Information: Client version 10.0.42.24719  
 Information: Client time 2016-08-02T11:19:55-04:00  
 Information: Eclipse Cone Dose Calculation (Version 10.0.28)  
 Information: Using DCF protocol version 0.3  
 Information: Beam data directory:  
 \\VARIANDB.ARIA.LOCAL\DCF\$client\BeamData\CDC\_10.0.28.2  
 Information: Calculated martes, 02 de agosto de 2016 11:25:21 a.m.

## Prescription

**Dose / Fraction:** 2200.000 cGy  
**Number of Fractions:** 1  
**Total Dose:** 2200.000 cGy  
**Treatment Percentage:** 70.0% of maximum  
**Plan Normalization Method:** No plan normalization  
**Plan Normalization Factor:** 1.000  
**Total Weight:** 10.740  
**Total Rotation:** 1800°  
**3D Dose Maximum (Dmax):** 3142.857 cGy  
**Repeat Factor:** 432.256 cGy

## Parameters

NOTE:

- All angles and field sizes are given using Varian IEC scale.
- Isocenter coordinates are given using the Planning coordinate system.
- Calibration Factors have been corrected to take the absolute dosimetry measurement geometry into account.
- Reference Dose = Repeat Factor × Weight Factor / Average TMR
- Monitor Units = Reference Dose / Calibration Factor

## Treatment Fields

Field ID	Isocenter				Cone [mm]	Couch Rotation	Gantry		Weight Factor	Cal.Fact. [cGy/MU]	Aver. TMR	Ref.Dose [cGy]	Monitor Units	MU/Deg [MU/°]	Aver.d [mm]
	ID	X [cm]	Y [cm]	Z [cm]			Start	Stop							
Campo 1	1	-2.19	-1.63	-2.20	14	10.0°	20.0°	120.0°	1.004	0.888	0.646	671.533	756 MU	7.56	106.3
Campo 2	1				12	30.0°	20.0°	120.0°	0.989	0.871	0.618	691.973	794 MU	7.94	113.2
Campo 3	1				12	50.0°	20.0°	120.0°	0.989	0.871	0.593	720.733	827 MU	8.27	121.0
Campo 4	1				12	350.0°	240.0°	340.0°	1.004	0.871	0.799	543.065	623 MU	6.23	62.4
Campo 5	1				12	330.0°	240.0°	340.0°	0.989	0.871	0.774	552.460	634 MU	6.34	68.8
Campo 6	1				14	310.0°	240.0°	340.0°	0.989	0.888	0.712	600.661	676 MU	6.76	86.5
Campo 7	2	-2.24	-1.63	-1.52	5	10.0°	20.0°	120.0°	0.330	0.715	0.619	230.718	323 MU	3.23	104.7
Campo 8	2				5	30.0°	20.0°	120.0°	0.330	0.715	0.594	240.127	336 MU	3.36	111.7
Campo 9	2				5	50.0°	20.0°	120.0°	0.330	0.715	0.571	250.126	350 MU	3.50	119.4
Campo 10	2				5	350.0°	240.0°	340.0°	0.330	0.715	0.777	183.591	257 MU	2.57	61.0
Campo 11	2				5	330.0°	240.0°	340.0°	0.330	0.715	0.754	189.403	265 MU	2.65	67.1
Campo 12	2				5	310.0°	240.0°	340.0°	0.330	0.715	0.686	208.078	291 MU	2.91	84.5
Campo 13	3	-2.93	-1.63	-1.80	5	10.0°	20.0°	120.0°	0.466	0.715	0.598	336.542	471 MU	4.71	111.1
Campo 14	3				5	30.0°	20.0°	120.0°	0.466	0.715	0.574	350.608	490 MU	4.90	118.1
Campo 15	3				5	50.0°	20.0°	120.0°	0.466	0.715	0.556	361.916	506 MU	5.06	123.9
Campo 16	3				5	350.0°	240.0°	340.0°	0.466	0.715	0.806	249.704	349 MU	3.49	54.3
Campo 17	3				5	330.0°	240.0°	340.0°	0.466	0.715	0.786	256.011	358 MU	3.58	59.0
Campo 18	3				5	310.0°	240.0°	340.0°	0.466	0.715	0.715	281.491	394 MU	3.94	76.6

Cal.Fact. = Calibration factor.

Aver.d = Average depth.